DD Bio CW

DE Gebrauchsanweisung, **EN** Instructions for Use, **FR** Notice d'utilisation, **IT** Istruzioni per l'uso, **ES** Instrucciones de uso, **NL** Gebruiksaanwijzing, **EL** Οδηγίες χρήσης για το, **PT** Instruções de utilização, **SV** Bruksanvisning, **ET** Kasutusjuhend, **DA** Brugsanvisning, **FI** Käyttöohje, **PL** Instrukcja używania, **CS** Rady pro použití, **HU** Használati útmutató, **RO** Instrucțiuni de utilizare, **SK** Návod na použitie, **BG** Инструкция за употреба, **HR** Upute za uporabu, **GA** Treoracha úsáide, **LV** lietošanas pamācība, **LT** Naudojimo instrukcija, **SL** Navodila za uporabo, **MT** Struzzjonijiet għall-użu

| | | | |
|-----------|-----------------|---------------------------|----|
| DE | DEUTSCH | Gebrauchsanweisung | 3 |
| EN | ENGLISH | Instructions for Use | 5 |
| FR | FRANÇAIS | Notice d'utilisation | 7 |
| IT | ITALIANO | Istruzioni per l'uso | 9 |
| ES | ESPAÑOL | Instrucciones de uso | 11 |
| NL | NEDERLANDS | Gebruiksaanwijzing | 13 |
| EL | ΕΛΛΗΝΙΚΑ | Οδηγίες χρήσης για το | 15 |
| PT | PORTUGUÊS | Instruções de utilização | 17 |
| SV | SVENSK | Bruksanvisning | 19 |
| ET | EESTI | Kasutusjuhend | 21 |
| DA | DANSK | Brugsanvisning | 23 |
| FI | SUOMI | Käyttöohje | 25 |
| PL | POLSKI | Instrukcja używania | 27 |
| CS | ČEŠTINA | Rady pro použití | 29 |
| HU | MAGYA | Használati útmutató | 31 |
| RO | ROMÂNĂ | Instrucțiuni de utilizare | 33 |
| SK | SLOVENSKY | Návod na použitie | 35 |
| BG | БЪЛГАРСКИ | Инструкция за употреба | 37 |
| HR | HRVATSKI | Upute za uporabu | 39 |
| GA | GAEILGE | Treoracha úsáide | 41 |
| LV | LATVIEŠU VALODA | lietošanas pamācība | 43 |
| LT | LIETUVIŲ KALBA | Naudojimo instrukcija | 45 |
| SL | SLOVENŠČINA | Navodila za uporabo | 47 |
| MT | MALTI | Struzzjonijiet għall-użu | 49 |

DD Bio CW

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DD Bio CW

Cobalt based dental alloy for ceramic, type 4*

*According to ISO 22674:2016

The product is intended for dental laboratories for the creation of products that must be rigid and resistant, single crowns, bridges and superstructures on implants.

The alloy DD Bio CW is produced in conformity with standards ISO 9693:2020 and ISO 22674:2016, it is highly corrosion resistant, it does not contain Nickel and it is completely free of toxic elements like beryllium, indium and gallium.

DD Bio CW is a product intended to be used by qualified technician.

Intended purpose

Fabrication of metal base structures for metal-ceramic crowns and FPDs.

Intended patient population

This alloy is intended to be used for oral rehabilitation in edentulous or partially edentulous patients who suffer from loss of hard dental tissue. There are no limitations related to age, sex, clinical condition or other conditions, other than those possibly determined case by case by the professional user.

Intended users

The product is intended for professionals of the dental sector, in particular: prosthodontic laboratories, dental technicians, dentists. These professional figures shall have a high degree of specialisation and are perfectly aware of the destination of the product and the correct method of use for the realization of a final product and elaborate different methods to reach the best shape of the final prosthesis. Contact and/or manipulation by other people is not envisaged.

Instruction for use

Modeling

Modeling with a minimum thickness of 0.3mm on simple structures and 0.5mm in complex cases or with patients with bruxism. Avoid elbows.

Manufacturing

Separate the milled structures from the disk with tungsten carbide burs or with separating discs, clean with ultrasonic and acetone or bioalcohol to eliminate any oil residues from the working process. Continue manufacturing process using tungsten carbide burs.

Stress relief heat treatment

As a first step we recommend a thermal relaxation treatment at 1030° for 7 minutes in vacuum (from 1 to 3 elements), while if we have more than 4 elements, we increase to 12 minutes. Once removed from the oven, the product must be sandblasted using 110 to 150 µm aluminum oxide as disposable material with a pressure of 3-4 bar. Steam and check that the surface has a homogeneous grey surface. After cleaning, the framework should not be touched anymore and should be held with clamp forceps only.

Re-use

The device can be reused until the available material is exhausted.

Ceramization

For synthetic ceramics it is possible to use Bonding. Apply the ceramic, following its manufacturer's instructions. For natural/feldspar ceramics, it is recommended to fire for 1 min under vacuum at 980°C. Do not use bonding and use directly the opaque diluted with distilled H₂O. In both cases, perform slow and longterm cooling. It is recommended to open the oven below 500°C. For detailed information, please refer to the instructions for use of the veneering material. Please only use veneering materials that have compatible CTE with DD Bio CW.

Soldering

We suggest using MAGNUM SALDATURA Co to weld.

Waste disposal

The processing scrap must be disposed of as special waste in accordance with the EC directives 2008/98/CEE on waste, and 94/62/CEE on packaging and packaging waste and in compliance with national legislation in force on the subject.

Safety instructions

- Cobalt-Chrome-based alloys can seldom cause dermatitis on sensitive subjects. A Patch-Test is thus advisable.
- Metal dusts and fumes are dangerous for health. Use an adequate exhaust fans during milling and polishing.
- The alloy contains Cobalt, element classified as CMR (CAS Number: 7440-48-4).
- Before prosthesis application verify if other metal elements are in patient's oral cavity. Coexistence of different metals can cause a "pile" effect.
- We identify every batch with a number. We recommend to write it down in patient's file to allow its complete traceability.
- This product does not need any special preservation precautions. MESA keeps batch traceability for 15 years.
- The methods of preserving the products made with DD Bio CW are under the responsibility of professional operators (dental technicians). This product must be stored, transported and used at a temperature between -30°C and +50°C to ensure the integrity of the packaging and the instructions contained therein.
- MESA recommends that the patient be made aware of the possibility for dental alloys to affect MRI results.
- The product is sold non-sterile.


Users shall report any major accident occurred in relation to the use of this medical device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

A summary of safety and clinical performance (SSCP) is available on Eudamed.

| Percentage chemical composition (m/m) | |
|---------------------------------------|----|
| Co | 66 |
| Cr | 27 |
| Mo | 6 |
| Others: Si, Mn | |

| Physical and mechanical properties | |
|--|---------------------------------------|
| Density | 8.4 g/cm ³ |
| Solidus-liquidus temperature | 1307-1417°C |
| Melting point | 1470°C |
| Yield load strength (Rp 0.2) | 395 MPa |
| Percentage elongation at fracture | 11% |
| Modulus of elasticity | 233 GPa |
| Vickers hardness | 255 HV10 |
| Colour | White |
| Thermal expansion coefficient 25-500°C | 14.3·10 ⁻⁶ K ⁻¹ |
| Thermal expansion coefficient 25-600°C | 14.5·10 ⁻⁶ K ⁻¹ |
| Maximum firing temperature | 980°C |
| Suggested ceramics | VITA VM13 |
| Cytotoxicity test according to ISO 10993-5 | Passed |

Symbol Description

| | |
|---|--|
|  | Manufacturer |
|  | Date of manufacture |
|  | Distributor |
|  | CE Marking |
|  | Batch number |
|  | Catalogue number |
|  | Non-sterile |
|  | Temperature limitation |
|  | Consult instructions for use www.dentaldirekt.de/en/IFU |
|  | Medical device |
|  | Product UDI Code |
|  | Contains hazardous substances |

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